



**Quality, reconditioned  
medical equipment that you can trust.**



 **Pocket Nurse®**  
Medical Supplies for Education and Simulation™

**(800) 225-1600 | [www.PocketNurse.com](http://www.PocketNurse.com)**



Defibrillators and AED



EKG / ECG



Exam Tables



Feeding Pumps,  
SSDs and ESUs



Infusion Therapy



Patient Monitoring



Patient Transport  
and Comfort



Ventilators



## The Pocket Nurse® Difference

We don't just wipe it down, spray some paint on it, and ship it. Our products undergo a 20-point inspection and are beautifully crafted to like-new condition with a one-year warranty!

## Biomed Inspection

Our factory-trained and certified technicians inspect, test, and return medical equipment to OEM standards.

## Cleaning, Disassembly, and Repair

We disassemble each product, thoroughly clean and disinfect, and replace worn or broken parts.

## Reconditioning

We professionally refinish our products in industrial paint booths using automotive-grade paints and spray guns for quality that you can see and feel.

## Two-Year Warranty



**stryker**



**Dräger**

**MAQUET**

**PHILIPS**

**WelchAllyn**

**Baxter**



**Hill-Rom**







## Zoll M Series Defibrillator

Lightweight, compact, and rugged defibrillator system. This defibrillator is half the weight and size of any competing unit in its class, and features tremendous value at an economical price. An 11 lb system that's packed with power. Easy to lift, position, and carry, it will never be left behind. Built-in AC power.

- Completely Reconditioned with One Year Parts & Labor Warranty
- Available in 3 lead or 12 lead EKG Configurations
- Available in Monophasic or BiPhasic Configurations
- Built-In recorder/printer for accurate documentation
- Field Emission Display for wider viewing angles than any LCD display. High contrast and visibility allow you to easily view display up to 60 degrees
- Compact size and light weight (11 pounds) allow for easy use and multiple applications
- Includes AC Power Cord / Battery Charger, rechargeable battery, Multi-Function Cable, Printer, one set of electrode pads, and one pad of printer paper
- Can be used with optional hard paddles (not included) in lieu of Multi-Function Cable and electrode pads



## Physio-Control, Inc. Lifepak 12 Defibrillator

Multi-parameter device that combines semi-automated and manual defibrillation with capnography, external pacing, 12-lead electrocardiography and other monitoring functions. Its platform is upgradeable, allowing customers to add new parameters as their needs change and as new features are introduced.



## Zoll E Series Defibrillator

Rugged and ready defibrillator intended for use by EMS personnel. ZOLL's Rectilinear Biphasic™ waveform (RBW) was designed specifically for external defibrillation to control for variations in patient impedance. Only the ZOLL RBW has demonstrated statistical clinical superiority to monophasic waveforms in peer reviewed randomized controlled trials.



## Physio-Control, Inc. Lifepak 20 Defibrillator

Compact and lightweight defibrillator that offers flexible and effective tools to meet diverse patient care needs across different settings. Using ADAPTIV™ biphasic technology, the LifePak 20 Defibrillator / Monitor offers both manual and AED modes for hospital and clinic settings. The system offers both monitoring and therapeutic functions.



## Zoll R Series Defibrillator

Truly two devices in one. To meet the needs of BLS-trained responders, it powers on as an AED. But with the push of a button, it converts to a full-featured manual defibrillator for ALS responders. Using standard wireless technology, it can automatically relay this information to be viewed on a PC or mobile device.



## ECG/EKG Machines

### MIPM TeslaDuo MRI Monitor

SKU: MIPM-TeslaDuo-Mon

Represents the latest standard in MRI pulse oximetry. TeslaDUO was specially designed for monitoring patients during MRI examinations. It meets the requirements of Radiologists and Anaesthesiologists without compromising image quality. The TeslaDUO continuously measures your patient's oxygen saturation and pulse as well as non-invasive blood pressure during MRI examinations. The TeslaDuo is FDA approved and safe for use as directed with MRI equipment.

- Brand new with 2 Year OEM Warranty
- Available in 3 different configurations: SPO2, NIBP or Combined with both SPO2 and NIBP
- 7" color touch screen display for good readability of the monitoring values, even from the control room
- Wireless SPO2 Sensor
- Infrared Remote Control
- USB Interface provides for easy data transfer



#### GE DASH 3000 Patient Monitor

8.4" Display. Easy-to-use interface and color-coded patient cables/connectors. Standard Monitor features NIBP, RESP, SPO2, Temp, and ECG, Optional CO2 and IBP monitoring available. Optional wireless LAN. Lithium-Ion internally or externally rechargeable battery(ies). 4-Channel Integrated Thermal Printer.



#### GE Procure B40 Patient Vitals Monitor

Provides clinical excellence in patient monitoring from GE Healthcare. Provide comprehensive, portable monitoring solutions at the bedside. Capable of monitoring adult, paediatric, and neonatal patients, these monitors are designed with advanced clinical measurement technologies, including proven GE algorithms.



#### GE DASH 4000 Patient Monitor

10.4" large-format, full-color active matrix display. Easy-to-use interface and color-coded patient cables/connectors. Optional wireless LAN. Long-lasting batteries power Dash 4000 for extended use.



#### GE MAC 1200 Diagnostic EKG

Compact, comprehensive solution with practical features that clinicians at smaller hospitals, clinics and physician offices will appreciate. Designed for acquisition with or without analysis, the system can be configured as a standard electrocardiograph without analysis for facilities with minimal ECG requirements.



#### GE DASH 5000 Patient Monitor

12.1" large-format, full-color active matrix display. Easy-to-use interface and color-coded patient cables/connectors. Optional wireless LAN. Long-lasting batteries. 3 waveform screen view. Optional CO2 module and extension rack. Both mainstream and low-flow sidestream CO2 measurement options.



#### GE MAC 5500 ECG / EKG

15 lead analysis with interpretation (12 SL EKG analysis program with Gender Specific Software). The information is displayed on a 10.4" backlit color AM-LCD display. Has the ability to provide full page printouts. With a single button for operating it can print, store and transmit data almost instantly. Optional barcode scanner.





## GE Carescape V100 Vital Signs Monitor with Rolling Stand

Low-acuity patient monitor that supports you in making fast, quality care decisions. The CARESCAPE V100 monitor is a reliable, easy-to-use vital signs monitor that helps you care for many patients efficiently. Its speed, accuracy and mobility allow you to collect the right information at the point of care, and present it wherever it is needed, so you can make fast, quality care decisions.

- Completely Reconditioned with One Year Parts & Labor Warranty
- Loaded Vital Signs Monitor – Measures BP, Temp, and SPO2
- Optional rolling stand, BP cuff/hose, SPO2 probe/cable, and Temperature probe
- Large, Bright LED Displays give you readings at a glance
- One step Start/Stop button
- Integrated Printer
- Battery automatically recharges when plugged in



## Datascorp Passport II Patient Monitor

Dedicated function keys Colorful waveforms with automatic display configuration Standard features include six waveforms, large numerics, 3 or 5-lead ECG, lead-selectable respiration, non-invasive blood pressure, temperature, and Masimo SpO2.



## GE Corometrics 259 / 259 CX Fetal Monitor

GE-Exclusive Smart BP. Minimizes monitor interactions by using our one touch hot keys. DINAMAP® SuperSTAT® non-invasive blood pressure technology Integrated Masimo SET or Nellcor maternal pulse oximetry. Selectable 3 Lead maternal ECG provides clinicians with QRS wave amplitude and meaningful diagnostic data.



## Welch Allyn 420 SPOT Vital Signs Monitor

Lightweight and portable. Big, bright, easy-to-read LCD screen. Measures NIBP, Pulse, Temp, and SPO2. SureTemp® Plus thermometry Optional Masimo® or Nellcor® pulse oximetry (SpO2). IR connectivity to interface with computerized patient records.



## Welch Allyn Spot LXI Vital Signs Monitor

SureBP measures blood pressure in about 15 seconds. SureTemp Plus Thermometer or Braun ThermoScan PRO 4000 Ear Thermometer. Masimo or Nellcor Pulse Oximetry. Bigger, brighter LCD display. Integrated long-life battery enhances portability and continuous operation.



## GE Corometrics 129 Fetal Monitor

Smart BP feature – delays auto determination during contractions. Vital signs history feature saves maternal information. Keyboard option for detailed annotations Interface with 340 Telemetry and Hydrotherapy when used with the Nautilus Watertight Transducers Interface capability with the OxiFirst® Fetal Pulse Oximetry system.





## Infusion Therapy

### Alaris Medley Infusion System

The BD Alaris™ Pump module is a large volume infusion pump that helps continuously or intermittently deliver fluids, medications, blood and blood products to adult, pediatric or neonatal patients. Clinicians can attach up to four IV infusion modules, allowing four independent infusions on a single PC unit.

Select your PC module, then select the infusion pump module that best meets your needs.



#### Alaris Medley 8015 PC Module Color Display

Modular smart IV medication safety system that can help protect every infusion to guard you and your patients from medication errors. The Alaris Medley Point of Care System combines a single point-of-care unit featuring Guardrails Suite MX safety software with IV / LVP, PCA Syringe, Syringe, respiratory monitoring (SpO2 and EtCO2), and barcoding capabilities.

#### Alaris Medley 8100 Infusion IV Pump

Drug Calculation allows entry of drug dose and flow rate. Downstream Occlusion Detection system. KVO Rate Adjust.  
Flow Rate Range: 0.1-99.9 ml/hr in 0.1 ml increments  
100-999 ml/hr in 1 ml increments  
Rate Accuracy:  $\pm 5\%$  at rates between 1 and 999 ml/hr  
 $\pm 5.5\%$  at rates less than 1 ml/hr

#### Alaris Medley 8110 Syringe Pump Module

Common user interface—same front panel keys as other modules  
Plunger drive head protected—extended handle to protect plunger.  
Flow Rate Range:  $\pm 2\%$  of full scale plunger  
100-999 ml/hr in 1 ml increments

#### Alaris Medley 8120 PCA Syringe Pump Module

Common user interface—same front panel keys as other modules  
Includes locking cabinet, keyset, and Patient Bolus Cord.  
Flow Rate Range:  $\pm 2\%$  of full scale plunger  
100-999 ml/hr in 1 ml increments

#### Alaris Medley 8300 EtCO2 Module

Displays trend real-time patient respiratory rate, EtCO2 and no breath alarm status. Uses Oridion Microstream technology with Filterline disposables (sold separately) to offer real time respiratory status.

#### Alaris Medley 8015 LS Next Gen PC Module Color Display

Provides a common user interface for programming infusions and continuous SpO2 / EtCO2 monitoring. This helps to reduce complexity at the point of care. Each Alaris Medley 8015 LS PC Module features a larger display and can support up to four Infusion Modules with a single power cord. The 8015 LS Module features a large, easy to read color display, and back-lit keypad for programming values. Programming screen displays drug or IV fluid name, dose or rate and volume infused. Real time battery gauge display lets clinicians know the amount of run time available for their current infusions. Alaris Medley 8015 LS Next Generation PC Modules are Guardrails compatible, and Educational customers have the option of having an Educational Training Dataset added for training purposes.





## Alaris Medsystem III Triple Channel Infusion Pump

Triple Channel Infusion System in a very small package. Onboard Rechargeable NiCad Battery.



## B Braun Infusomat Space Infusion Pump

Independent Modularity, Easily Stackable, KeyGuard™: Non-numeric keypad, Extensive Drug Library, Wireless Integration.



## B Braun Outlook Series Infusion Pumps

Real-time display of dose limit values during infusion for titration safety. Active alphabetic matrix.



## B Braun Outlook Series Infusion Pumps

Real time display of dose limit values during infusion for titration safety. KeyGuard™: Non-numeric keypad.



## B Braun Perfusor Space Infusion Pump

Alarms are audible and color-coded for easy recognition. Syringes automatically recognized for fast set-up.



## B Braun Vista Basic Infusion Pump

DEHP-free sets protect patients. Needle-free sets protect both the patient and the clinician.



## Baxter / Sigma Spectrum Infusion Pump

Bright, large-type, color display. Baxter's standard Clear-Link and Clear-Flo IV administration sets. Opt Wireless.



## Baxter 6201 Flo-Gard Infusion Pump

Single Channel Volumetric I.V. Pump. Front Panel Lock-out. Flow Check Occlusion Alarm.



## Baxter 6301 Flo-Gard Dual Channel Infusion Pump

Dual Channel Volumetric I.V. Pump. Flow Check Occlusion Alarm. Programmed delivery profiles.



## Baxter AS50 Syringe Pump

Flow Rate Accuracy:  $\pm 3\%$  of full scale plunger travel  
Volume Accuracy:  $\pm 3\%$ , or 0.007" of travel.



## Baxter InfusOR Syringe Pump

Interchangeable Smart Labels. Large Easy to Read Dials. Battery life is typically 150 hours.



## Covidien / Kendall Kangaroo 924 Enteral Feeding Pump

Six different audible and visual function alarm. Highly visible LED display. (1 to 300 mL/hour) and VTBD (1 to 2000 mL).



## Covidien / Kendall Kangaroo ePump Enteral Feeding

Feed Rate – 1-400 ml/hr in 1 ml increments  
Flush Rate – 10-500 in 1 ml increments.



## Covidien / Kendall Kangaroo Joey Ambulatory Enteral Feeding Pump

Pump accuracy is  $\pm 7\%$  High IPX Rating. DEHP-Free Pump Sets. Anti-Free Flow Pump Sets.



## Hospira Lifecare PCA 3 Syringe Pump

Wireless connectivity. Locking Medication Cabinet. Power Cord, Pole Clamp, Keypad, and Patient Bolus Cord Included.



## Hospira Lifecare PCA Syringe Pump

Barcode Reader. Locking Medication Cabinet. Power Cord, Pole Clamp, Keypad, and Patient Bolus Cord Included.



## Hospira Plum A Plus Infusion Pump

Single Channel or Triple Channel. Integral set-based, free-flow protection. PlumSet™ Administration Set with CLAVE®.



## ICU Medical / Hospira Plum 360 Infusion Pump

PlumSet™ cassette that features air management. EMR Ready. Automated programming and documentation. Wireless.



## Medfusion 2010 i Syringe Pump

Flow rate: 0.01-378.0 ml/hr dependent upon syringe size  
Bolus rate: up to 378.0 ml/hr dependent upon syringe size.



## Medfusion 3500 Syringe Pump

PharmGuard® Medication Software. FlowSentry™ monitor. Customizable drug alerts.



## Moog Curlin 4000 CMS Ambulatory IV Pump

SKU: MOO-CURL-4KC

Flow Rates 0.1mL to 400mL/hr. Real-time Line Pressure Display. Volume Limits 1.0mL to 9,999mL.



## Moog Curlin 6000 CMS Ambulatory IV Pump

Real-time Line Pressure Display. Highly Visible Large Fonts. Titration in continuous and PCA Modes.



## Patient Transport and Comfort

### Hill-Rom Durastar P8005

The ideal choice for the ER, OR, PACU, Recovery, or for outpatient settings. This rugged stretcher offers outstanding mobility, comfort, and excellent quality. Accommodates virtually all patient conditions, as well as your needs in this uncertain environment.

- Completely reconditioned with one year parts & labor warranty
- **Best value and most reliability in class**
- Standard 700 lb capacity
- Optional features include Intellidrive, onboard scale, push handles, and IV pole
- Available in 26" and 30" widths
- Retractable side rails with full extension bumper
- Hilow and Trendelenburg/reverse control pedals and Brake/steer pedals – head and foot
- Low height for safer stretcher ingress and egress – best in class



### Hill-Rom Transtar P8000 Stretcher

**Best value and most reliability in class**  
500 lb capacity and optional 700 lb capacity  
Available in 26" and 30" widths. Retractable side rails with full extension bumper. Hilow and Trendelenburg/reverse control pedals and Brake/steer pedals.



### Hill Rom 8050 OB / GYN Stretcher

**Best value and most reliability in class**  
700 lbs weight capacity. Sliding litter enables quick and easy patient positioning. Next generation Steering Plus – system for crisp cornering and control. Dual-locking, four wheels brakes for strong, stable platform.



### Pedigo 540 Stretcher

500-pound patient weight capacity. Precision casters that lock from any corner. Pneumatic assist positioning. Hands-free hydraulic height adjustment. Trendelenburg and reverse Trendelenburg. Tuck-away side rails for gapless positioning.



### Stryker 5050 / 5051 Stretcher Chair

Pneumatic backrest maneuvers from Zero to 90 degrees. Dual hydraulic jacks and four wheel, steel-ring brakes. Foot section can be moved dependently or independently of backrest.



### Stryker Prime Series Stretchers

700 lb weight capacity. Integrated Utility Tray and thermo-formed base with extra capacity storage tray  
4 Wheel, center-locking brake system with dual end activators. Dual, side mounted foot control hydraulics with uni-lower pedal. Dual end brake / steer controls.



### Ferno 35A Lightweight Ambulance Cot

Lightest cot available with a 500 pound (227 kg) load capacity. Foot-end lead handle and side release handle  
Two loading heights. One-hand swing-down sidearms. Includes patient pad and straps.



### Ferno PROFlexx 35X Ambulance Cot

All-Metal Frame with Built-in Shock Frame. 5-Position Drop Frame with One Handed drop feature. One Handed Safety Bail Release. Telescoping Side Lift Handles. Wider Bed Surface and Head Frame  
Foot-End Lead Handle.



### Stryker MX Pro Ambulance Cot

Heavy-duty oversize rear wheels. X-frame design. Tough, highly visible yellow powder coating. MX-PRO's standard one-hand release. Fold down side rails. Weight Capacity: 500 – 650 lbs, depending on model.





## Hill-Rom P1600 Advanta

Easy-grip siderail controls that are conveniently located in a central location where both patient and caregiver interaction naturally occurs. Four corner brake and steer system allows for easy mobility, and a hard pan mattress surface coupled with a CPR quick release function is available to allow for quick CPR positioning along with a removable headboard to assist in CPR positioning. The P1600 Advanta Bed also features ergonomically upgraded head and footboards, as well as siderails. This makes the bed much more comfortable to work on and around for patients and care providers.

- Completely reconditioned with one year parts & labor warranty
- Includes Brand new Hospital-Grade Mattress
- **Our Most Popular / Best Value Hospital Bed**
- Sleep surface: 36" x 83"
- Overall Length: 91" Retracted Length: 78.5"
- Automatic Contour Up (Maximum): 15"
- 3-Year Warranty



### Hill Rom 1100 Series Advance Bed

**Monet Medical recommended - best low-cost value!**  
Includes Brand new Hospital-Grade Mattress  
Total Length 91" Retracted Length 78.5"  
Mattress Dimensions: 83" X 35" X 6"  
Maximum Width (siderails up): 42.75"



### Hill-Rom Care Assist ES Bed

SlideGuard® bed articulation technology with autocontour. Integrated scale. Point-of-Care® nurse controls. One-button dining chair®. Optional IntelliDrive® powered transport or central 5th steering wheel.



### Hill-Rom Affinity I / II OB/GYN Birthing Bed

Manual Operation in case of power failure. Built-in foot supports adjust vertically and horizontally to position feet. Automatic night light increases safety for patients and nurses. Hinged seat and back sections keep mattress from bunching up.



### Hill-Rom P3200 VersaCare Bed

Bed goes into full chair position. Easy-to-grip side rails. Fixed headboard stays the same - the optimal height for pushing. Hands-free foot controls to adjust bed. Optional 3-Level Bed Exit Alarm.



### Hill-Rom Affinity III OB/GYN Birthing Bed

Tuck away side rails for minimum transfer gap. Includes new mattress covers, calf supports, and new Placenta Basin. Manual Operation and Battery Back Up in case of power failure. Built-in foot supports adjust vertically and horizontally to position feet.



### Stryker Secure II Hospital Bed

Dual Pedestal design. CPR Release. In-rail Fowler/Knee Gatch control. Electronic function lock-out control. Four-wheel, steel-ring brake system with centrally-located activation pedals. Backrest angle indicator. Locking-caster wheel steering with 4 traction sockets.



### Hill-Rom Affinity IV OB/GYN Birthing Bed

Tuck away side rails for minimum transfer gap. Manual Operation and Battery Back Up System in case of power failure. Built-in foot supports adjust vertically and horizontally to position feet. Hinged seat and back sections keep mattress from bunching up.



### Midmark Pedigo 500 Pediatric Crib / Stretcher

The unique rounded corners reduce the overall size of the 500, yet provide the same amount of useful area on the 26 1/2" x 55 3/4" mattress. The two-post support system provides 97% access to the patient by the medical team.



## Exam Tables, Power Chairs, and OR Tables

### Midmark Ritter 204 Exam Table

Features the latest in design, comfort and functionality in a manual exam table, it truly offers efficient patient care. Midmark found that you need an exam table that allows easy access for patients and makes them feel secure. The Midmark 204 Exam Table also features smart options like three side drawers for optimal supply management, modern styling to coordinate with your office, an updated, ergonomic, comfortable top design, and additional weight capacity to handle a wider range of patients.

- Completely Reconditioned with One Year Parts & Labor Warranty
- Brand new Upholstery – over 50 different available colors
- One-Piece Pass-Through Seamless Drawers
- Stylish and Easy to Clean
- Removable Top System requires no tools
- Enlarged Footstep
- Storage Access and Volume
- Increased Weight Capacity



### Midmark Ritter 75 Table/Chair

Power procedure table/chair. It features adjustable positioning for: OB/GYN, colposcopy, general procedures, ENT, proctology, urology, laser surgery, etc. It features controls which automatically and separately adjust and position the height, tilt, back section, and leg section. There are easy manual adjustments for head rest and proctological positioning. Importantly, there is an automatic return positioning button.



### Midmark 104 Exam Table

Provides exceptional value in a manual exam table, and it truly offers efficient patient care. Monet Medical completely reconditions these tables and offers your choice of over 50 different colors of upholstery. This Examination Table offers easy patient positioning, spring-assisted backrest, hide-away stirrups which pull out and lock in numerous positions, non-skid locking footstep, and five large drawers for storage.



### Skytron 6500 O.R. Table

The exclusive 180° rotating top allows unobstructed C-arm and microscope coverage. Additionally, it eliminates conflicts between surgeon and anesthesiologist during head-end procedures, and provides ample leg room for surgeons who prefer to sit during procedures.



### Steris AMSCO 4085 O.R. Table

The 4085 general surgical table meets your needs with a versatile, high-productivity general surgical table specifically designed to provide outstanding value. Now you can get a table that offers excellent imaging and patient posturing capabilities all in one.



### Steris / AMSCO 3085 SP O.R. Table

One of the most reliable and popular surgical tables. Trusted for years by healthcare facilities across the world, the 3085 is a cost-effective solution for the general OR. The AMSCO 3085 SP surgical table deftly handles virtually all the surgical specialties you expect a general surgical table to handle,





## Carefusion / VIASYS AVEA Ventilator

Provides comprehensive critical care ventilation for a wide range of patients with complex respiratory diseases. In the ever-changing world of healthcare, you face more challenging patients and more pressure to deliver outcomes. The AVEA was designed with you in mind. Its comprehensive array of clinical tools and advanced monitoring allow you to guide your patient through the continuum of care simply and efficiently with a focus on improving outcomes.

- Completely Reconditioned with One Year Parts & Labor Warranty
- Integrated Rolling Cart/Stand
- Humidifier: Active and passive humidification
- Leak compensation
- Mode type: A/C, SIMV, CPAP/PSV, NPPV, nasal CPAP/IMV
- Breath type: APRV/BiPhasic, Volume, Pressure, TCPL, PRVC, Volume Guarantee
- Apnea backup: Volume, pressure



## Drager Evita XL Ventilator

Integrated Trolley / Stand Design makes for easier movement throughout care facilities. Mask ventilation with NIV plus available in all modes with comprehensive monitoring and alarm adaptation. Screen: 15" TFT color touch-screen.



## Phillips Respironics Trilogy 100 / 200 / 202 Portable Ventilator

Compact design with long-life internal battery for intra-hospital transport. Supports active and passive circuits for invasive and noninvasive ventilation. Innovative leak compensated volume control ventilation.



## Puritan Bennett PB840 Ventilator

Integrated Rolling Stand/Cart. Standard 9.4" Graphic Display Unit or Optional 10.4" Graphic Display Unit available. Optional on-board air compressor for use in facilities without house air. Standard Adult / Pediatric Ventilator with many available options including Neo-Natal.



## Respironics Esprit Ventilator

Includes Rolling Stand/Cart. Available with optional Patient Breathing Circuit and Test Lung. Pressure Support, Pressure Control Ventilation. Fully integrated microprocessor control for system control and monitoring. Pressure Control and Non-Invasive Ventilation Settings. Graphical User Interface (GUI) Touchscreen.



## Siemens Servo - I Ventilator

Optional Compressor Available. Optional Patient Breathing Circuit and Test Lung available. Integrated Rolling Stand. Invasive or Non-Invasive Intubation. Non-Invasive Intubation has leakage compensation. Transportation capable.



## Respironics BiPap Vision Ventilator

Includes Rolling Stand and support arms. Integrated display screen displays real-time graphics. CPAP and S/T modes. Oxygen control allows 21-100% concentrations from without attaching auxiliary oxygen to the mask.





Please feel free to call anytime.  
To view our full line of products, visit us online.

**(800) 225-1600 | [www.PocketNurse.com](http://www.PocketNurse.com)**

Availability of equipment and services shown in this booklet are subject to change without notice.  
Please contact your Account Manager to assist you in finding the solutions that best meet your needs.